ENGLISH

STAAR® nanoPOINT™ Single Use Lens Injector System
For use with STAAR® AFINITY™ COLLAMER® SINGLE PIECE IOL


DEVICE DESCRIPTION
The sterile, single use nanoPOINT™ INJECTOR SYSTEM is a device to implant STAAR AFINITY™ Collamer Single Piece IOL (Collamer IP IOL). The delivery system provides a tubular pathway through an incision over the iris, to the center of the pupil. The IOL may then be gently injected from the tube into the eye. The nanoPOINT™ INJECTOR SYSTEM consists of three parts:

1. The nanoPOINT™ INJECTOR (Fig. 1) with silicone lip.
 2. The loading unit (Fig. 2).
 3. The nanoPOINT™ cartridge (Fig. 3).
- INDICATION FOR USE**
The nanoPOINT™ INJECTOR SYSTEM is indicated for loading and injecting the Collamer IP IOL into the capsular bag following extracapsular cataract extraction.
- CONTRAINDICATIONS**
Use of the STAAR Injector with other intraocular implants, or with cartridges that have not been tested and approved by STAAR Surgical Company.
- INSTRUCTIONS FOR USE**
1. Open the injector system packaging in the sterile area and review the sterile contents.
2. The nanoPOINT™ cartridges are equipped with guiding rails. This allows rapid and precise loading of the lens. Hold the loading unit with the nanoPOINT™ cartridge in such a manner that the rear cartridge flange can be guided with the

index finger of the left hand. Position the rear cartridge flange so that the distance between the guiding rails corresponds approximately to the size of the lens eye (Fig. 4).
3. Push the lens onto the guide rails from the rear and position it in the middle of the loading chamber (Fig. 5).
4. Apply viscoelastic material directly to the opening of the cartridge barrel, as well as to the loading chamber below the lens (Fig. 6).
5. Close the cartridge on the loading unit (the lens automatically inserts itself into the cartridge) until the locking mechanism engages. Carefully remove the cartridge from the loading unit (Fig. 7a, 7b).
CAUTION: The lens should be injected within two minutes after loading. Viscoelastic material may be gently lubricated if allowed to start too long while exposed to air.
6. Inject the cartridge into the opening at the front end of the injector and push it to the very front position (Fig. 8a). Carefully remove the plunger forward and guide the silicone lip

into the loading chamber of the cartridge. Push the lens out of the chamber into the tapered lip of the cartridge until the leading plate hinge reaches the lip of the cartridge (Fig. 8b).
CAUTION: Retract plunger and then push forward again. The lens enters the eye as the lens is always grasped correctly.
GUARANTEE AND LIABILITY LIMITATION
The manufacturer guarantees that this product has been manufactured with the appropriate care and assumes no liability for side effects or resulting damages, losses or costs that may occur as the direct or indirect use of this product. Manufacturer's liability is restricted to the performance of repairs resulting from product defects which are clearly not the result of incorrect handling.
EXPIRATION DATE
The expiration date on the nanoPOINT™ INJECTOR SYSTEM packaging is the sterility expiration date. The sterility is assured if the tray seals are not punctured or damaged prior to the expiration date. The injector system cannot be used past the indicated expiration date.

REPORTING
Adverse events and / or potentially life-threatening complications that may reasonably be regarded as implant and injector related and that were not previously expected in nature, severity or degree of incidence should be reported to STAAR Surgical at:
USA Phone: (800) 352-7842
Fax: (800) 952-4923
HOW SUPPLIED
The nanoPOINT™ INJECTOR SYSTEM is supplied in sterile sealed trays. The packaging is sterilized by an appropriate method and should be opened only under sterile conditions.
RETURN POLICY FOR nanoPOINT INJECTOR SYSTEM
Contact STAAR Surgical.
STORAGE
Store at room / ambient temperature.

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
CAUTION: Federal (USA) law restricts the sale of this device to, or upon the order of, a physician.

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
Directions for use
Gebrauchsanweisung
Istruzioni d'uso
Mode de empleo

Manufactured by:
MEDICEL AG
Luchten 1262
CH-9427 Wolfhalden
SWITZERLAND

Tel. +41 71 727 10 50
Fax +41 71 727 10 55
Info www.medicelag.com
www.medicelag.com

Distributed by:
STAAR Surgical
1911 Walker Avenue
Monrovia, CA 91016
USA
Tel: (800) 352-7842
Fax: (800) 952-4923

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